

Refraction Station Specifications:

Frame Material	Steel
Wood Top	1.5-inch-thick solid maple
Weight -- Frame and Top	85 lbs.
Tabletop Dimension	48" X 31"
Sliding Instrument Table Size	13" X 14"
Refractor arm height	24.65" – 44.95"
Vertical Movement	Dual motors
Lifting capacity	250lbs
Lifting speed	1.1" /sec
Noise level	45 dB
Height range	27.45" – 47.75"
Outlet voltage	110 – 240V
Anti – collision	Yes
Electrical outlets	3
Programmable	up to 4 height settings
Refractor arm height	24.65" – 44.95"
Included accessories	1 Refractor bar, 1 cable channel, 5mm wrench